MAGFIT IP-F Series




MAGFIT IP-F series is compatible with Frialit2 Implant system

- •This is compatible with 3.8mm, 4.5mm type
- •There are three length to be able to adapt to each patient's gingiva condition.
- •It has TiN-coating at the attractive face not to be worn
- It is tightened by using the JIS standard 1.3mm hexagonal driver.

Specifications

Product Lineup



Keeper Construction (mm)	0.7 0.7 1.7 03.8/04.5	0.7 3.0 3.8/04.5	0.7 5.0 03.8/04.5
Assembly[Dome]	17type	30type	50type
DXD	FDD17S FDD17L	FDD30S FDD30L	FDD50S FDD50L
Attractive force	600 gf		

Materials

- ·Keeper Screw : AUM20
- $\cdot \text{Keeper}$ attractive face is coated by TiN for wear resistance.
- ·Abutment ring : Ti

Compatible implant models

Frialit2(Zive) ϕ 4.5, ϕ 3.8

Tightening torque Recommendation

20 to 25N·cm

Tools

To tighten the keeper, it requires a specialized tool or JIS Standard 1.3mm hexagonal driver.