

MAGFIT™ IP-I Series



MAGFIT IP-I series is compatible with ITI Implant system.

- There are three length to be able to adapt to each patient's gingival condition.
A taper is added to the fitting part of the fixture to prevent loosening of the screw
- It has TiN-coating at the attractive face not to be worn.
- It is tightened by using the JIS standard 1.3mm hexagonal driver.

Specifications

Product Lineup

Keeper Construction	Magnetic Assembly [Flat]			Keeper Construction	Magnetic Assembly [Dome]		
	14 type	30 type	45 type		14 type	30 type	45 type
Magnetic Assembly [Flat]	14 type	30 type	45 type	Magnetic Assembly [Dome]	14 type	30 type	45 type
DXFL	IFD14	IFD30	IFD45	DXD	IDD14	IDD30	IDD45
Attractive force	750 gf			Attractive force	600 gf		
SXFL	IFX14	IFX30	IFX45				
Attractive force	550 gf						

Materials

- Keeper Screw: AUM20
- Keeper attractive face is coated by TiN for wear resistance.
- Abutment ring: Ti

Compatible implant models

- Branemark: Standard Implant, No.043.030S-035S, 043.131S-135S, 043.230S-234S(SLA)
Standard plus implant, No.043.050S-054S, 043.151S-154S, 043.250S-234S(SLA)
 - Swiss Plus: ϕ 4.8mm D platform, No.SP8-14, SPWB8-14, OPB8-14, OPWB8-14
- For more detail information, please make an enquiry

Tightening torque Recommendation

25 to 30N·cm

Tools

To tighten the keeper, it requires a specialized tool or JIS Standard 1.3mm hexagonal driver.