

MAGFIT™ IP-I Series






MAGFIT IP-I series is compatible with ITI Implant system.

- There are three length to be able to adapt to each patient's gingival condition.
A taper is added to the fitting part of the fixture to prevent loosening of the screw
- It has TiN-coating at the attractive face not to be worn.
- It is tightened by using the JIS standard 1.3mm hexagonal driver.

Specifications

Product Lineup

Keeper Construction	Magnetic Assembly [Flat]			Keeper Construction	Magnetic Assembly [Dome]		
	14 type	30 type	45 type		14 type	30 type	45 type
	IFD14	IFD30	IFD45		IDD14	IDD30	IDD45
Attractive force	750 gf			Attractive force	600 gf		
	IFX14	IFX30	IFX45				
Attractive force	550 gf						

Materials

- Keeper Screw: AUM20
- Keeper attractive face is coated by TiN for wear resistance.
- Abutment ring: Ti

Compatible implant models

- Branemark: Standard Implant, No.043.030S-035S, 043.131S-135S, 043.230S-234S(SLA)
Standard plus implant, No.043.050S-054S, 043.151S-154S, 043.250S-234S(SLA)
 - Swiss Plus: ϕ 4.8mm D platform, No.SP8-14, SPWB8-14, OPB8-14, OPWB8-14
- For more detail information, please make an enquiry

Tightening torque Recommendation

25 to 30N·cm

Tools

To tighten the keeper, it requires a specialized tool or JIS Standard 1.3mm hexagonal driver.