# MAGFIT IP-V Series



MAGFIT IP-V series is compatible with Replace select Implant system.

- •This is compatible with only regular platform.
- •There are three length to be able to adapt to each patient's gingiva
- •It has TiN-coating at the attractive face not to be worn.
- •It is tightened by using the JIS standard 1.3mm hexagonal driver.

## **Specifications**

#### **Product Lineup**

Keeper Construction	0.5	0.5	0.5
Magnetic Assembly [Flat]	30 type	40 type	55 type
DXFL	VFD30	VFD40	VFD55
Attractive force	750 gf		
SXFL	VFX30	VFX40	VFX55
Attractive force		550 gf	

Keeper Construction	04.7	0.7	0.7
Magnetic Assembly [Dome]	30 type	40 type	55 type
7	VDD30	VDD40	VDD55
DXD			
Attractive force	600 gf		

#### Materials

- ·Keeper Screw: AUM20
- · Keeper Screw attractive face is coated by TiN for wear resistance.
- ·Abutment ring: Ti

#### Compatible implant models

NB Replace Select Straight / Tapered Regular Platform (Ø4.3mm) For more detail information, please make an enquiry.

### Tightening torque Recommendation

25 to 30N·cm

#### Tools

To tighten the keeper, it requires a specialized tool or JIS Standard 1.3mm hexagonal driver.